

Intracellular Staining (IFN Gamma, IL-4, etc) Protocol

DO NOT USE THIS PROTOCOL FOR FOXP3 STAINING

Example:

	FL1	FL2	FL3	FL4	FL5
12	Thy1.2	isotype		CD8	CD4
13	Thy1.2	IFN Gamma		CD8	CD4
14	Thy1.2	IL-4		CD8	CD4

1. Equally aliquot cells into FACS tubes in preparation for staining
2. Take cells that are to be stained intracellularly, combine into one FACS tube, and centrifuge. Stain other cells using normal protocol
3. Make incubation media for intracellular tubes while cells are spinning (200 μ L needed per FACS tube)
 - a. Using 1640 complete media, add 0.1 μ L / 1mL media ionomycin (Sigma # I0634) and 0.0065 μ L / 1mL media PMA (Sigma # P8139; for example, I add 0.32 μ L ionomycin and 0.021 μ L PMA to 3.2mL 1640 CM*)
 - b. Briefly vortex incubation media
4. Once cells spin down, discard supernatant and resuspend in 200 μ L incubation media per FACS tube (using above example, resuspend in 600 μ L incubation media)
5. Pipet 200 μ L into separate wells of a round-bottom (U bottom) 96-well plate (using above example, 3 wells are used)
6. Place plate in cell incubator for one hour
7. Add 1 μ L of GolgiPlug (BD # 555029) per 1mL media to each well (in example, 0.2 μ L of GolgiPlug is added to each individual well since each contains 200 μ L media)
8. Put plate back in incubator for 2.5 hours
9. After time expires, mechanically dislodge cells via pipetting, place cells back in separate FACS tubes, spin, and stain as usual

*CM = Complete Media = RPMI 1640 w/ 10% FBS, 1% Pen/Strep (complete RPMI 1640)